



## Habib™ Hexablate for ablation of primary and secondary cancers of the liver



- Comprises ergonomic handle, aspiration lines and electrode array
- Outer array allows for the creation of a cylindrical sealed area
- Central needle allows for both aspiration and sealing of inferior and superior areas
- Provides treatment for liver tumours unsuitable for resection



- Bipolar
- Aspirator channel
- Switch box allows sequential switching between outer array and central needle
- Ablation zone created at periphery of tumour
- Removal of vapour and aspirate reduces possibility of contamination
- 10mm and 20mm open and laparoscopic versions

- CE marked
- Cleared for marketing by FDA
- Supplied sterile
- Ready to use
- Single use product

**EMcision Limited**  
**13-15 Wallace Way**  
**Hitchin**  
**Herts**  
**SG4 0SE**  
**UK**

**Website: [www.emcision.com](http://www.emcision.com)**  
**Phone: +44 (0) 1462 423594**  
**Contact: [info@emcision.com](mailto:info@emcision.com)**